Low Pressure Connector Tubing with Prime Tube Tubulure de raccord basse pression avec tubulure de purge Tubo conector de baja presión con tubo de cebado

Instructions for Use / Mode d'emploi / Instrucciones de Uso

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© 2017-2019 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Stellant, and MEDRAD Stellant are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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A glossary of the symbols used on the above listed sterile disposable products can be found on page 2 of these Instructions for Use.

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ENGLISH

Introduction: Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner.

Indications for Use: The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single-use on one patient only with MEDRAD® Stellant injectors.

Contraindications: These devices are not intended for multiple patient use, drug infusions, chemotherapy, or any other use for which the device is not indicated.

Restricted Sale: Rx Only

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

Warnings

This product is intended for single use only. Do not re-sterilize, reprocess, or reuse. The disposable devices have been designed and validated for single-use only.

Re-use of the single-use disposable devices poses risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned and re-sterilized.

Do not use if sterile package is opened or damaged.

Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.


Patient or operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures of 400 psi (2750 kPa). Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures.

Contamination may occur if the end of the luer connector is touched. Do not touch the end of the luer connector.

Cautions

Refer to the injector Operation Manual and syringe Instructions for Use for further instructions.

Remove Dust Caps


1. Remove the dust cap on the luer fittings.

Cautions

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.

Attach Tubing

1. Ensure all air is purged from the syringe.
2. Attach the connector tube to the syringe, 1/4 turn to 1/2 turn maximum.

NOTE: If using a connector tube with T-connector, attach the straight portion of the T-connector to the contrast (Syringe A) and the extension to the saline (Syringe B). If the T-connector is connected to Syringe B, priming will not completely fill the connector tubing.

Attach Prime Tube


1. Connect the Prime Tube to the patient end of the tubing set.

Warnings

Clear all air from the syringes and tubing before connecting a patient to the injector. Air embolization can cause death or serious injury to the patient.

Automatic prime may not function properly and an air injection may result if SJS-LP-60-T-J is used without additional manual priming. Use manual priming to ensure all air is removed from the tubing.

Prime Tubing & Check for Air

1. Press **Prime** on the injector.
2. Ensure all air is expelled.
3. Rotate the injector head downward.

Remove & Discard Prime Tube


1. Remove and discard the Prime Tube.

Warnings

Reuse of this product may result in biological contamination, product degradation, and/or product performance issues. Properly discard disposable items after use, or if there is any possibility that contamination may have occurred.

Connect to Patient

1. Press **I Checked for Air**, confirming that the operator checked the fluid path for the presence of air.
2. Connect the tubing to the patient.

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